

FUMEHOOD/
LEMARI ASAM

MYCO7-FH120

Efisien, elegan, dan aman untuk semua kebutuhan laboratorium.

Aplikasi Ideal

- Laboratorium pendidikan dan penelitian
- Industri farmasi, kosmetik, dan kimia
- QC/QA laboratorium
- Laboratorium rumah sakit dan universitas



Fume hood ekonomis dengan perlindungan maksimal. Kombinasi material plywood + HPL + PP yang tangguh, anti-korosi, dan elegan untuk kebutuhan laboratorium modern.

Keunggulan:

- ✓ Anti-korosi & mudah dirawat
- ✓ Desain modern & ekonomis
- ✓ Operasi senyap, efisien, dan stabil
- ✓ Komponen modular, mudah diservis
- ✓ Produk lokal, standar global

- Garansi 1 tahun
- Sertifikat uji fungsi manufaktur
- Manual book




FUMEHOOD/LEMARI ASAM

MYCO7-FH120



SPECIFICATIONS :

Body dimensions	: 120 x 70 x 220 cm (W.D.H)
External material	: HPL - laminated plywood
Internal material	: polypropylene (PP) sheet
Working plate	: granite / phenolic resin
Front sash door	: 6 mm clear glass with manual slide up system
Fan	: 180 watt centrifugal blower PP polypropylene 450 Pa, 1100 CMH
Noise level	: < 75 dB (low noise level)
Exhaust air velocity	: > 0.4 m/s
Exhaust duct	: PVC / flexible pipe
Working lamp	: 28 watt LED downlight
Switches	: LED, fan, & power switch
Electrical socket	: 2 unit waterproof socket
Accessories	: sink, water tap



FUMEHOOD AIRFLOW SYSTEM

MYCO7

PRECISION • SAFETY • INNOVATION