

FUMEHOOD/
LEMARI ASAM

MYCO7-FH150

Fume hood berukuran 150 cm dengan desain luas dan elegan. Terbuat dari plywood berlapis HPL (luar) dan polypropylene (PP) (dalam) yang tahan terhadap bahan kimia korosif.

Aplikasi Ideal

- Laboratorium pendidikan dan penelitian
- Industri farmasi, kosmetik, dan kimia
- QC/QA laboratorium
- Laboratorium rumah sakit dan universitas



DIMENSI :
150 x 80 x 220 cm (P . L . T)

Perlindungan Optimal. Desain Lebih Luas. Kinerja Lebih Stabil.



Keunggulan:

- ✓ Ruang kerja lebih luas (150 cm)
- ✓ Material kuat & tahan bahan kimia
- ✓ Sistem ventilasi aman
- ✓ Lampu LED hemat energi
- ✓ Lengkap dengan sink & water tap
- ✓ Desain modular, mudah dirawat

- Garansi 1 tahun
- Sertifikat uji fungsi manufaktur
- Manual book




FUMEHOOD/LEMARI ASAM MYCO7-FH150



SPECIFICATIONS :

- Body dimensions : 150 x 80 x 220 cm (W.D.H)
- External material : HPL - laminated plywood
- Internal material : polypropylene (PP) sheet
- Working plate : granite / phenolic resin
- Front sash door : 8 mm clear glass with manual slide up system
- Fan : 180 watt centrifugal blower PP polypropylene 450 Pa, 1100 CMH
- Noise level : < 75 dB (low noise level)
- Exhaust air velocity : > 0.4 m/s
- Exhaust duct : PVC / flexible pipe
- Working lamp : 28 watt LED downlight
- Switches : LED, fan, & power switch
- Electrical socket : 2 unit waterproof socket
- Accessories : sink, water tap



FUMEHOOD AIRFLOW SYSTEM

MYCO7

PRECISION • SAFETY • INNOVATION