

INKUBATOR CO2 MYCO7-IC1200

- Solusi Laboratorium Modern dengan Kapasitas 120 Liter dan Kontrol Presisi
- Dirancang untuk Stabilitas Suhu, CO2 Akurat, dan Lingkungan Kultur Steril
- Produk Lokal Pertama untuk Kultur Sel dan Penelitian Biologi



DIMENSI P.L.T :
600 X 650 X 850 mm

Stabil. Akurat. Steril.

Fitur Unggulan:

- Produk lokal berkualitas internasional
- Suhu & CO2 stabil, hasil riset lebih akurat
- Filter HEPA untuk lingkungan kultur steril
- Layanan purna jual cepat & mudah
- Cocok untuk laboratorium riset, universitas, hingga industri bioteknologi

- Garansi 1 tahun
- Sertifikat uji fungsi manufaktur
- Manual book



INKUBATOR CO2 MYCO7-IC1200



SPECIFICATIONS



External Material : Steel plate coating
Internal Material : Mirror stainless steel
Shelves : Adjustable, made of stainless steel
Capacity : 120 Liters
External Dimensions : 600 mm (W) x 650 mm (D) x 850 mm (H)
Internal Dimensions : 500 mm (W) x 500 mm (D) x 550 mm (H)
Temperature Range : RT +5°C to 60°C
Temperature Accuracy : ± 0.5°C



Heating system : Direct heating
CO2 Concentration : 0 - 5%
Range
CO2 Accuracy : ± 0.2%
CO2 Sensor : Infrared (IR Sensor) for precise control
Humidity Range : 70 - 90% RH



Alarm System : Over temperature, sensore failure,
and door open alerts

Display : LCD display
Electrical Power : 750W
Electrical Voltage : 220V /50Hz
Filter : HEPA filter to maintain a sterile
environment



Additional Accessories : - Adjustable shelves (3-5 levels)
- Lockable wheels for easy mobility

Gross Weight : ± 80 kg

Applications : Ideal for cell culture, biological research,
and microbiology laboratories

MYCO7

PRECISION • SAFETY • INNOVATION

