

MINI FUMEHOOD /LEMARI ASAM

MYCO7-FH80

memberikan perlindungan optimal terhadap uap asam, gas beracun, maupun debu kimia berbahaya.

Aplikasi Ideal

- Laboratorium pendidikan dan penelitian
- Industri farmasi, kosmetik, dan kimia
- QC/QA laboratorium
- Laboratorium rumah sakit dan universitas



Mini Fume Hood MYCO7-FH80 dirancang khusus untuk laboratorium dengan ruang terbatas yang tetap membutuhkan sistem ventilasi aman dan efisien.

Keunggulan:

- Portabel & hemat ruang
- Material anti-korosi
- Suara rendah, nyaman digunakan
- Pencahayaan LED hemat energi
- Soket waterproof
- Kaca Tempered Slide-Up

- Garansi 1 tahun
- Sertifikat uji fungsi manufaktur
- Manual book




MINI FUMEHOOD / LEMARI ASAM MYCO7-FH80



SPECIFICATIONS :

| | |
|----------------------|---|
| Body dimensions | : 80 x 70 x 130 cm (W.D.H) |
| External material | : HPL - laminated plywood |
| Internal material | : polypropylene (PP) sheet |
| Working plate | : granite / phenolic resin |
| Front sash door | : 5 mm clear glass with manual slide up system |
| Fan | : 80W polypropylene Fan 8", 2450 RPM, 800 CMH, 220V |
| Noise level | : < 75 dB (low noise level) |
| Exhaust air velocity | : > 0.4 m/s |
| Exhaust duct | : PVC / flexible pipe |
| Working lamp | : 18 watt LED downlight |
| Switches | : LED, fan, & power switch |
| Electrical socket | : 2 unit waterproof socket |



FUMEHOOD AIRFLOW SYSTEM

MYCO7

PRECISION • SAFETY • INNOVATION