

# MINI LAMINAR AIR FLOW MYCO7-LM60

Ciptakan ruang kerja steril dan bebas kontaminasi dengan LAF MYCO7-LM60. Dirancang untuk efisiensi penggunaan ruang tanpa mengurangi performa sterilisasi.

## Aplikasi Utama:

- Kultur jaringan tanaman & jamur
- Inokulasi bakteri & mikroba
- Peracikan obat farmasi
- Sterilisasi alat & bahan laboratorium
- Operasi alat elektronik sensitif



**DIMENSI :**  
60 x 60 x 180 cm ( P . L . T )

## Steril. Efisien. Presisi.

### Kenapa Memilih LM60?

Ukuran kompak, cocok untuk laboratorium kecil

Material berkualitas, tahan lama, dan mudah dibersihkan

Harga efisien & terjangkau

- Garansi 1 tahun
- Sertifikat uji fungsi manufaktur
- Manual book




# MINI LAMINAR AIR FLOW MYCO7-LM60



## SPECIFICATIONS :

- Laminar type : horizontal
- Body dimensions : 60 X 60 X 110 cm (W.D.H)
- Full dimensions : 60 X 60 X 180 cm (W.D.H)
- Material : HPL - laminated plywood
- Working plate : stainless steel 304
- Front sash door : 5 mm thickness clear/tempered glass with manual slide system
- Filter : HEPA filter H13 efficiency of 99,95% @ MPPS (0,1–0,3  $\mu$ m)
- Fan : 65 watt blower RQA 20060 HBL
- Noise level : < 75 dB (low noise level)
- Air flow velocity : > 0.3 m/s
- Sterilization lamp : 10 watt UVC germicidal lamp 254 nm
- Work lamp : 9 watt LED lamp
- Switches : LED, UV, fan & power switch



Air Circulation and Filtration Pattern  
in Horizontal Laminar Air Flow

# MYCO7

PRECISION • SAFETY • INNOVATION