

1. Pilih Jenis Aliran Udara

VERTIKAL – Udara bersih turun dari atas ke bawah.

➔ Cocok untuk: kultur jaringan, mikrobiologi, farmasi.

✓ Lebih aman dari risiko kontaminasi dari operator.

HORIZONTAL – Udara bersih mengalir dari belakang ke depan.

➔ Cocok untuk: elektronika presisi, optik, pengujian alat.

✓ Aliran udara lebih stabil dan lembut.

2. Tentukan Ukuran & Material Sesuai Aktivitasmu

Model	Ukuran Kabinet (mm)	Material	Ideal Untuk
LM60	600 × 600 × 1100	Plywood + HPL	Praktikum siswa, lab sekolah, pelatihan dasar
LM80	800 × 700 × 1300	Plywood + HPL	Kultur jaringan dasar, riset pendidikan
LM120	1200 × 750 × 1400	Plywood + HPL	Laboratorium kampus & pelatihan riset
LM600	600 × 650 × 1100	Stainless steel 304 + powder-coated steel	Praktek siswa / lab kecil
LM800	800 × 700 × 1300	Stainless steel 304 + powder-coated steel	Kultur jaringan dasar
LM1200	1200 × 750 × 1400	Stainless steel 304 + powder-coated steel	Riset & rutinitas laboratorium
LM1500	1500 × 850 × 1400	Stainless steel 304 + powder-coated steel	Dua operator ringan
LM2000	2000 × 900 × 1400	Stainless steel 304 + powder-coated steel	Produksi besar / multi operator
S1200	1200 x 750 x 1400	Full stainless steel 304	Riset & rutinitas laboratorium
S2000	2000 x 900 x 1400	Full stainless steel 304	Produksi besar / multi operator

MYCO7 Laminar Air Flow Product Index (2026–2027)

HALAMAN

MODEL

UKURAN (mm)

MATERIAL

04 - 05	LM60	600×600×1100	Plywood + HPL
06 - 07	LM80	800×700×1300	Plywood + HPL
08 - 09	LM120	1200×750×1400	Plywood + HPL
10 - 11	LM600	600×650×1100	Stainless steel 304 + powder-coated steel
12 - 13	LM800	800×700×1300	Stainless steel 304 + powder-coated steel
14 - 15	LM1200	1200×750×1400	Stainless steel 304 + powder-coated steel
16 - 17	LM1500	1500×850×1400	Stainless steel 304 + powder-coated steel
18 - 19	LM2000	2000×900×1400	Stainless steel 304 + powder-coated steel
20 - 21	S1200	1200×750×1400	Full stainless steel 304
22 - 23	S2000	2000×900×1400	Full stainless steel 304



Note: seluruh unit Laminar Air Flow MYCO7 dapat dikustomisasi sesuai kebutuhan pelanggan — termasuk tipe aliran udara, ukuran kabinet, dan desainnya.

MYCO7 Laminar Air Flow Price List (2026–2027)

MODEL	UKURAN (mm)	MATERIAL	HARGA
LM60	600x600x1100	Plywood + HPL	39.700.000
LM80	800x700x1300	Plywood + HPL	48.700.000
LM120	1200x750x1400	Plywood + HPL	69.700.000
LM600	600x650x1100	Stainless steel 304 + powder-coated steel	74.700.000
LM800	800x700x1300	Stainless steel 304 + powder-coated steel	96.700.000
LM1200	1200x750x1400	Stainless steel 304 + powder-coated steel	132.700.000 (touchscreen)
LM1500	1500x850x1400	Stainless steel 304 + powder-coated steel	165.700.000 (touchscreen)
LM2000	2000x900x1400	Stainless steel 304 + powder-coated steel	223.700.000 (touchscreen)
S1200	1200x750x1400	Full stainless steel 304	154.700.000
S2000	2000x900x1400	Full stainless steel 304	243.700.000

Note. Harga untuk unit dengan spesifikasi atau ukuran custom akan disesuaikan berdasarkan desain dan kebutuhan klien.

